

HEMODIALYSIS LABORATORY TEST REQUEST

ACCOUNT NAME AND ADDRESS

10329 Stony Run Lane Ashland, VA 23005

Phone (804) 365-3000

Toll Free (800) 888-8061

Fax (804) 365-3002

DATE SHIPPED

MACHINE NUMBER (SAMPLE ID)

FACILITY NAME

SAMPLE DATE

TYPE OF WATER (CLASSIFICATION) **DO NOT SEND DIALYSATE**

SPECIAL INSTRUCTIONS AND/OR UNUSUAL CONDITIONS:

CHECK ONE:

☐

AAMI ANALYSIS

ANSI/AAMI RD62-2006

Aluminum Chromium Potassium

Antimony Copper Selenium

Arsenic Fluoride Silver

Barium Lead Sodium

Beryllium Magnesium Sulfate

Cadmium Mercury Thallium

Calcium Nitrate Zinc

☐

LAL/Endotoxin (Special Kit)

☐

HPC

☐

TOC

☐

Other

CONTACT NAME

TELEPHONE NUMBER

EMAIL/ FAX NUMBER

HEMODIALYSIS WATER SAMPLE KIT INSTRUCTIONS

Each kit contains one bottle one zip lock plastic bag, and one pre-labeled mailer.

1. Run water tap or access line for two (2) minutes.
2. Rinse bottle with a small amount of water to be sampled.
3. Fill bottle completely with water sample, cap tightly, and place filled bottle in the plastic bag.
4. Fill out the Hemodialysis Laboratory Test Request form.
5. Place bottle enclosed in plastic bag and completed test request inside mailing container and mail.

*Ship LAL/HPC analysis cold but not frozen overnight, first morning delivery to preserve the integrity of the sample.

If you have any questions, please call Customer Services at 800-888-8061.

PLEASE COMPLETE A SEPARATE FORM FOR EACH WATER SAMPLE