

Order Number	Company		Contact Name	
	Job No #		Email Address	
ATO808	Phone Number		Date required	
	Lab Pickup		Delivery Address	

1L PASSIVATED CANISTER

TELESCOPIC CANISTER		
Quantity	Canister Volume	Flow Rate Sampling Interval (Flow Rates)
	1L	

6L PASSIVATED CANISTER

Quantity	Canister Volume	Flow Rate Sampling Interval (Flow Rates) / Certification
	6L	

PASSIVE SAMPLERS - RAD145®

Quantity	Item
	RAD145® - Thermal Desorption, Yellow Diffusion Body, Support Plate
	Outdoor Shelter

ACTIVE SAMPLERS - T017

Quantity	Item
	TO-17 Thermal Desorption Tube

SUMMA CANISTER – ACCESSORIES

Quantity	Item
	Duplicate T-piece
	2-Propanol (40 mL)
	1/4" Swagelok nuts + Ferrules
	Quick Connects
	3-Way Valve
	Additional Pressure Gauge
	Triplicate t-piece

Equipment List (LABORATORY ONLY)

Certification Number	Canister ID	Flow Controller ID	Integrity Check
Additional Requests			