



Your complete partner for radon measurement services

From our Global Centre of Excellence for radon expertise, you can order both passive and electronic radon detectors.



• Long-term measurement:

Indoor radon testing for periods ranging from 60 days to one year.

• Short-term measurement:

Indoor radon testing conducted over a period of 7 to 14 days.

• Soil Radon Measurement

:From 5 to 7 days



radon.eurofins.se



info.radon.se@etn.eurofins.com



010-490 84 80



• **Ecosense EcoTracker:**

A sniffer device that quickly and easily identifies radon variations during mitigation work.



• **Ecosense RadonEye Pro:**

Logger designed for follow-up radon measurements in workplaces and schools.



radon.eurofins.se



info.radon.se@etn.eurofins.com



010-490 84 80