

Eurofins Milk Testing Denmark A/S
Ladelundvej 85
DK-6600 Vejen
Tlf. +45 70 22 42 86

Calibration Samples



Production week: 23 - 2024
Production date: 03-06-2024

Reference: Combifoss, Reference Manual Calibration

Uhomogenized raw milk		
Cell 1	Cell 2	Cell 3
207.000	378.000	609.000

Date: 06-06-2024

A handwritten signature in black ink, appearing to read "Bente Østergaard", written over a horizontal line.

Version 1

CRM Approver, Bente Østergaard

Contact: calibration@ftn.eurofins.com