



eurofins

Allergen Test Codes

Eurofins Food Testing UK



Allergen	Test code (food and rinse waters)	ELISA test kit	Lower LOQ (mg/kg)	LOQ max (mg/kg)	Test Code (swabs)	Detection limits (ng/ml) swabs
Wolverhampton laboratory						
Almond	UD0FT	SENSISpec® ELISA Almond	0.4	10	UD0FU	20-500
BLG	UD0F5	SENSISpec ELISA Beta-Lactoglobulin	0.1	4	UD0F6	0.5-20
Casein	UD0F2	SENSISpec ELISA Casein	0.2	6	UD0F4	10-300
Casein (meat, sausage, ham)	UD0F3	SENSISpec ELISA Casein	1	30	UD0F4	10-300
Milk proteins	UD0HR	SENSISpec ELISA Milk	0.4	10	UD0HP	20-500
Crustacean (tropomyosin)	UD0G1	SENSISpec ELISA Crustaceans (Tropomyosin)	0.020 (20µg)	0.400 (400µg)	UD0G2	1-20
Whole egg protein	UD0FH	Egg (Ovalbumin) ELISA kit II (Morinaga BioScience, Inc.)	0.31	20	UD0FI	0.78-50
Gluten	UD0EY (for gluten free products only)	SENSISpec INgezim Gluten R5	2	50	UD0EZ	4-100
Gluten	UD20A (for samples containing gluten)	SENSISpec INgezim Gluten R5	2	50	UD0EZ	4-100
Hazelnut	UD0FV	SENSISpec ELISA Hazelnut	1	40	UD0FW	50-2000
Hazelnut	UDA09	RIDASCREEN® FAST Hazelnut (R-Biopharm®)	2.5	20	UDA11	125-1000
Mustard	UD0FZ	SENSISpec ELISA Mustard	2	60	UD0G0	100-3000
Peanut	UD0FR	SENSISpec ELISA Peanut	1	40	UD0FS	50-2000
Sesame	UD0EF	SENSISpec ELISA Sesame	2	30	UD0EG	100-1500
Soya	UD0G3	RIDASCREEN FAST Soya (R-Biopharm)	2.5	20	UD0G4	25-200
Walnut	UD0FX	SENSISpec ELISA Walnut	2	60	UD0FY	100-3000

SENSISpec® is a brand of Gold Standard Diagnostics.

Page 1 of 2



eurofins

Allergen Test Codes



Allergen	Test code (food and food products)	ELISA test KIT	Lower LOQ (mg/kg)	LOQ max (mg/kg)	Test code (swabs)	Detection limits (ppm/swab) swabs
Linlithgow laboratory						
Gluten	ZB01H	Ridascreen® Gliadin (R-Biopharm®)	5	80	ZB01I	5 - 80
Soya	ZB01M *	Ridascreen Fast Soya (R-Biopharm)	2.5	20	ZB01P	2.5 - 20
Whole egg protein	ZB026 * (and ZB01K)	Ridascreen Egg (R-Biopharm)	0.12	0.98	ZB027 (and ZB01L)	0.12 - 0.98

*Accredited to food and food products (excluding alcoholic and soft drinks); swabs are accredited.

Blue: Internal test codes required for analytical testing but not reported in Analytical Reports.

Red: Swabs units are ppm/swab; different from our Wolverhampton laboratory.

Allergen	Test code (food and food products and environmental swabs)	PCR test kit	LOD (mg/kg)
Heathrow laboratory			
Soya**	CC011 + CC010	SureFood® allergen soya (R-Biopharm®)	1
Mustard**	CC012 + CC010	SureFood allergen mustard (R-Biopharm)	1
Celery**	CC013 + CC010	SureFood allergen celery (R-Biopharm)	1.5

** Test codes required include CC010 for DNA Extraction Advanced using SureFood® PREP Advanced (R-Biopharm®)

Eurofins Food Testing UK Ltd offers allergen testing and consultancy services to support accurate food labelling and effective allergen risk management for food businesses operators.

Using validated laboratory methods, we can help verify allergen presence, investigate potential cross-contact, and support due diligence against regulatory and customer requirements, giving you confidence in your allergen controls and labelling decisions.



FTUK148FY1