

Accreditation



The Deutsche Akkreditierungsstelle attests with this **Accreditation Certificate** that

Eurofins Food Integrity Control Services GmbH
Berliner Straße 2, 27721 Ritterhude

operates a testing laboratory that fulfills the requirements according to DIN EN ISO/IEC 17025:2018 for those conformity assessment activities specified in detail in the annexes listed below. This includes additional existing legal and normative requirements for the testing laboratory including those in relevant sectoral schemes, provided that these are explicitly confirmed in the annexes listed below.

D-PL-20965-01-01 Valid from: 19.01.2026

The management system requirements of DIN EN ISO/IEC 17025 are written in the language relevant to the operations of testing laboratories and they conform to the principles of DIN EN ISO 9001.

This accreditation was issued in accordance with Art. 5 Para. 1 Sentence 2 of Regulation (EC) 765/2008, after an accreditation procedure was carried out in compliance with the minimum requirements of DIN EN ISO/IEC 17011 and on the basis of a review and decision of the appointed accreditation committees.

This accreditation certificate only applies in connection with the notice of 19.01.2026. It consists of this cover sheet, the reverse side of the cover sheet and the corresponding annex

Registration number of the accreditation certificate: **D-PL-20965-01-00**

Berlin, 19.01.2026 Barbara Tyralla | Head of Technical Unit

Translation issued: 06.02.2026

This accreditation certificate was issued by the Deutsche Akkreditierungsstelle GmbH (DAkkS). It is digital sealed and valid without signature. It reflects the status as indicated by the date of issue. The current status of any valid and surveyed accreditation can be found in the directory of accredited bodies maintained by Deutsche Akkreditierungsstelle GmbH (www.dakks.de).

This document is a translation. The definitive version is the original German accreditation certificate.

See notes overleaf

Deutsche Akkreditierungsstelle GmbH

Office Berlin
Spittelmarkt 10
10117 Berlin

The Deutsche Akkreditierungsstelle GmbH (DAkkS) is the entrusted national accreditation body of the Federal Republic of Germany according to § 8 section 1 AkkStelleG in conjunction with § 1 section 1 AkkStelleGBV. DAkkS is designated as the national accreditation authority by Germany according to Art. 4 Para. 4 of Regulation (EC) 765/2008 and clause 4.7 of DIN EN ISO/IEC 17000.

Pursuant to Art. 11 section 2 of Regulation (EC) 765/2008, the accreditation certificate shall be recognised as equivalent by the national authorities within the scope of this Regulation as well as by the WTO member states that have committed themselves in bilateral or multilateral mutual agreements to recognise the certificates of accreditation bodies that are members of ILAC or IAF as equivalent.

DAkkS is a signatory to the multilateral agreements for mutual recognition of the European co-operation for Accreditation (EA), International Accreditation Forum (IAF) and International Laboratory Accreditation Co-operation (ILAC).

The up-to-date state of membership can be retrieved from the following websites:

EA: www.european-accreditation.org

ILAC: www.ilac.org

IAF: www.iaf.nu

Deutsche Akkreditierungsstelle

Annex to the accreditation certificate D-PL-20965-01-01 according to DIN EN ISO/IEC 17025:2018

Valid from: 19.01.2026

Date of issue: 19.01.2026

This annex is part of the accreditation certificate D-PL-20965-01-00.

Holder of the accreditation certificate:

**Eurofins Food Integrity Control Services GmbH
Berliner Straße 2, 27721 Ritterhude, Germany**

with the location

**Eurofins Food Integrity Control Services GmbH
Berliner Straße 2, 27721 Ritterhude, Germany**

The testing laboratory meets the requirements of DIN EN ISO/IEC 17025:2018 to carry out the conformity assessment activities listed in this annex. The testing laboratory meets additional legal and normative requirements, if applicable, including those in relevant sectoral schemes, provided that these are explicitly confirmed below.

The management system requirements of DIN EN ISO/IEC 17025 are written in the language relevant to the operations of testing laboratories and they conform to the principles of DIN EN ISO 9001.

*This annex to the certificate was issued by Deutsche Akkreditierungsstelle GmbH (DAkkS) and is digitally sealed.
This annex to the certificate is only valid together with the written accreditation certificate and reflects the status as indicated by the date of issue. The current status of any valid and surveyed accreditation can be found in the directory of accredited bodies maintained by Deutsche Akkreditierungsstelle GmbH (www.dakks.de).*

Annex to the accreditation certificate D-PL-20965-01-01

Tests in the fields:

Consumer health protection
Flexible scope of accreditation:

Within the test areas indicated, the testing laboratory is permitted, without being required to inform and obtain prior approval from DAkkS

[Flex A] To use standardised or equivalent test methods listed here with different issue dates.

[Flex C] To modify, develop or further develop test methods.

The testing laboratory has an up-to-date list of all test methods within the flexible scope of accreditation. The list is publicly available on the website of the testing laboratory.

In the case of flexibility in accordance with categories B and C, all test procedures, including their issue statuses, within the flexible scope of accreditation and, if applicable, their modification/restriction, are only shown in the published list of test procedures within the flexible scope of accreditation. The same applies to issue statuses in the case of flexibility in accordance with category A.

Consumer health protection
Physical, physico-chemical and chemical analysis
Foodstuffs

Type of test Test method (detector)	Analyte/ measurand	Matrix / test item	Short title Standard/normative or in-house method	Modification and/or restriction	[Flex]
Chromatography Liquid chromatography with conventional detectors (UVD, ELSD, FLD, RID)	Ingredients	Honey, bee products, plant- based foodstuffs			C
Chromatography Liquid chromatography with mass selective detectors (MS/MS, HR- MS)	Ingredients, additives, pesticide residues, residues, contaminants	Honey, bee products, plant- based foodstuffs			C
Chromatography Gas chromatography with conventional detectors (FID)	Residues	Honey			C
Chromatography Gas chromatography with mass selective detectors (MS)	Pesticide residues, residues	Honey			C

Valid from: 19.01.2026

Date of issue: 19.01.2026

Page 2 of 4

This document is a translation. The definitive version is the original German annex to the accreditation certificate.

Annex to the accreditation certificate D-PL-20965-01-01

Type of test Test method (detector)	Analyte/ measurand	Matrix / test item	Short title Standard/normative or in-house method	Modification and/or restriction	[Flex]
Chromatography Ion chromatography with conventional detectors (CD, EC)	Ingredients	Foodstuffs			C
Chromatography Ion chromatography with mass- selective detectors (MS/MS)	Pesticide residues, contaminants	Bee products			C
Gravimetry	Water-insoluble substances	Honey	DIN 10743		A
Photometry	Ingredients	Honey			C
Refractometry	Water content	Honey	DIN 10752-2	Brix measurement, digital refractometer with automatic temperature compensation	A
Spectrometry Elemental analysis using isotope- ratio mass spectrometry (EA-IRMS)	Ingredients	Honey			C
Titrimetry	Conductivity, pH value, free acids	Honey	ICS SOP 510-06: 2021-04		-

Sensory analysis
Foodstuffs

Type of test Test method (detector)	Analyte / measurand	Matrix / test item	Short title Standard/normative or in-house method	Modification and/or restriction	[Flex]
Simple descriptive tests	Appearance, consistency, smell, taste	Honey	ICS SOP 520-02: 2024-01		-

Valid from: 19.01.2026

Date of issue: 19.01.2026

Page 3 of 4

This document is a translation. The definitive version is the original German annex to the accreditation certificate.

Annex to the accreditation certificate D-PL-20965-01-01

Visual analysis

Foodstuffs

Type of test Test method (detector)	Analyte / measurand	Matrix / test item	Short title Standard/normative or in-house method	Modification and/or restriction	[Flex]
Optical microscopy	Ingredients, impurities	Honey, plant- based foodstuffs			C

Abbreviations used:

DIN	Deutsches Institut für Normung e. V. – German Institute for Standardization
EN	European standard
ICS SOP XXX-XX	In-house method of Eurofins Food Integrity Control Services GmbH
IEC	International Electrotechnical Commission
ISO	International Organization for Standardization

Valid from: 19.01.2026

Date of issue: 19.01.2026

Page 4 of 4

This document is a translation. The definitive version is the original German annex to the accreditation certificate.